

Anigen Rapid FCoV Ab

General Information

Purpose : Detection of Feline Corona virus antibody

Principle : Immunochromatographic assay

Specimen : Whole blood, plasma or serum

Materials provided :

- Anigen Rapid FCoV Ab test Device
- Assay diluent buffer
- Anticoagulant tube
- Disposable capillary tube

Special Features :

- One-step testing procedure : Quickly & Accurately
- High sensitivity and specificity
- Sensitivity : 96.0 % vs. IFA
- Specificity : 97.9 % vs. IFA




How to do FIP diagnosis with Rapid FCoV Ab test

-Many healthy cats and the cats with conditions other than FIP are FCoV Ab positive.

-Some cats with effusive FIP appear to have low titres or are FCoV Ab negative.

-FCoV Ab testing should only be performed to diagnose FIP in conjunction with a compatible history, clinical signs, and examination of effusions or blood for high globulins and a low A:G ratio.

Test Procedure



Ordering Information

Cat. No.	Description	Type	Pack size
RC12-04	FIV Ab/FeLV Ag	Dual Device	10 Tests/Kit