

FPV Ag

General Information

Purpose : Detection of Feline Parvo virus antigen

Principle : Immunochromatographic assay

Specimen : Feline feces

Materials provided :


- Anigen Rapid FPV test Device
- Specimen tubes containing assay diluent buffer
- Sample collection swabs
- Disposable droppers

Special Features :

- One-step testing procedure : Quickly & Accurately
- High sensitivity and specificity
- Sensitivity : 97% vs. RT-PCR
- Specificity : 98.5% vs. RT-PCR



Test Procedure



Ordering Information

Cat. No.	Description	Type	Pack size
RG12-03	FPV Ag	Device	10 Tests/Kit