

Giardia Ag

General Information

Purpose : Detection of vegetative stages and cysts of Giardiasis

Principle : Immunochromatographic assay

Specimen : Canine feces

Materials provided :


- Anigen Rapid Giardia test Device
- Specimen tubes containing assay diluent buffer
- Sample collection swabs
- Disposable droppers

Special Features :

- Able to detect the pathogen in both dogs and cats
- One-step testing procedure : Quickly & Accurately
- High sensitivity and specificity
- Detection limit : 125 cysts/100ml of feces



Test Procedure



Ordering Information

Cat. No.	Description	Type	Pack size
RG18-04	Giardia Ag	Device	10 Tests/Kit