

UltraSlim® LED Illuminator



SLB-01W



SLB-01

Specifications

210 x 210 x 40 (DxWxH)

2.3Kg

470nm, Blue LED light source

Sensitivity: up to 0.5 ng

< 10 % CV (coefficient of variance)

Viewing surface(mm): 120D x 70W

Feature

Compact: Enlarge your laboratorial space


Light: Convenient to be moved around

Safe: No damage to you and your samples

Accurate: Highly detectable ability

uniform: Observe anywhere you need

User Friendly: Easy to view and cut gel



Much Cheaper, and Much Smaller

