

OrTube | Centrifugal Tubes



**Low Protein
Binding PP**

**Individually
Leak-Tested**



**Special rack hold-
ing either 20 x 50
ml tubes or 30 x 15
ml tubes, or a mix**

**Featuring Alpha-
numerical identifi-
cation**



**Engraved graduation on
the cone, allowing easy
reading of small volumes**

CATALOGUE #	DESCRIPTION	DIAMETER (MM)	MATERIAL	VOLUME (ML)	UNITS / BAG	UNITS / CASE
5540100	50ml Tube conical	28,3	PP	50	20	440
5540150	50ml Tube conical, non sterile	28,3	PP	50	20	440
5540200	50ml Tube skirt	28,3	PP	50	40	400
5540300	15ml Tube conical	16,1	PP	15	40	560
5540350	15ml Tube conical, non sterile	16,1	PP	15	40	560
5540400	15ml Tube round	16,1	PP	15	40	560
5040500	15ml Tube conical PS	16,1	PS	15	40	560
5540620	Rack				1	25

Other TC Products Available




OrFlask Tissue Culture Flask

Unique square and large opening.
Patented OrCap squeeze system.




OrDish TC Round Dishes

4 Portusions allow easy and safe handling.



OrPlate TC Test Plates

Deep inter well spacing to the outer edge preventing cross contamination.



OrScrapers TC Scrapers

Special blade design. Easy, safe Grip.
Reaches easily all flask surface.